

PR23A rabbit pAb

Cat No.:ES19106

For research use only

Overview

Product Name	PR23A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PR23A AA range: 198-248
Specificity	This antibody detects endogenous levels of human PR23A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃. Avoid repeated freeze-thaw cycles.
Protein Name	Proline-rich protein 23A
Gene Name	PRR23A
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	729627
Human Swiss-Prot Number	A6NEV1
Alternative Names	
Background	

