

# RAG1 rabbit pAb

Cat No.:ES19110

For research use only

## Overview

Product Name	RAG1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human RAG1 AA range: 894-944
Specificity	This antibody detects endogenous levels of human RAG1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	V(D)J recombination-activating protein 1 (RAG-1) (RING finger protein 74) [Includes: Endonuclease RAG1 (EC 3.1.-.-); E3 ubiquitin-protein ligase RAG1 (EC 6.3.2.-.)]
Gene Name	RAG1 RNF74
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	5896
Human Swiss-Prot Number	P15918
Alternative Names	
Background	

