



# G3ST1 rabbit pAb

Cat No.:ES19112

For research use only

## Overview

<b>Product Name</b>	G3ST1 rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB
<b>Species Cross-Reactivity</b>	Human;Mouse
<b>Recommended dilutions</b>	WB 1:1000-2000
<b>Immunogen</b>	Synthesized peptide derived from human G3ST1 AA range: 280-330
<b>Specificity</b>	This antibody detects endogenous levels of human G3ST1
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C . Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	Galactosylceramide sulfotransferase (GalCer sulfotransferase) (EC 2.8.2.11) (3'-phosphoadenosine-5'-phosphosulfate:GalCer sulfotransferase) (3'-phosphoadenylylsulfate:galactosylceramide 3'-sulfotransf
<b>Gene Name</b>	GAL3ST1 CST
<b>Cellular localization</b>	Golgi apparatus membrane ; Single-pass type II membrane protein .
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	
<b>Human Gene ID</b>	9514
<b>Human Swiss-Prot Number</b>	Q99999
<b>Alternative Names</b>	
<b>Background</b>	

