



GEMI4 rabbit pAb

Cat No.:ES19113

For research use only

Overview

Product Name	GEMI4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human GEMI4 AA range: 754-804
Specificity	This antibody detects endogenous levels of human GEMI4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Gem-associated protein 4 (Gemin-4) (Component of gems 4) (p97)
Gene Name	GEMIN4
Cellular localization	Cytoplasm. Nucleus. Nucleus, nucleolus. Nucleus, gem. Localized in subnuclear structures next to coiled bodies, called gems, which are highly enriched in spliceosomal snRNPs and in the nucleolus.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	50628
Human Swiss-Prot Number	P57678
Alternative Names	
Background	

