



MKRN3 rabbit pAb

Cat No.:ES19118

For research use only

Overview

Product Name	MKRN3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MKRN3 AA range: 237-287
Specificity	This antibody detects endogenous levels of human MKRN3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Probable E3 ubiquitin-protein ligase makorin-3 (EC 6.3.2.-) (RING finger protein 63) (Zinc finger protein 127)
Gene Name	MKRN3 D15S9 RNF63 ZNF127
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	7681
Human Swiss-Prot Number	Q13064
Alternative Names	
Background	

