

TALAN rabbit pAb

Cat No.:ES19119

For research use only

Overview

Product Name	TALAN rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TALAN AA range: 41-91
Specificity	This antibody detects endogenous levels of human TALAN
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Talanin
Gene Name	ZNF365 KIAA0844
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	22891
Human Swiss-Prot Number	Q70YC4
Alternative Names	
Background	

