

## CCD82 rabbit pAb

## Cat No.:ES19126

For research use only

## Overview

Product Name	CCD82 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
<b>Recommended dilutions</b>	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CCD82 AA range: 322-372
Specificity	This antibody detects endogenous levels of human CCD82
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20 $^\circ\!\mathrm{C}$ . Avoid repeated freeze-thaw cycles.
Protein Name	Coiled-coil domain-containing protein 82
Gene Name	CCDC82 HT025
<b>Cellular localization</b>	
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	79780
Human Swiss-Prot Number	Q8N4S0
Alternative Names	
Background	



ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C