



ANS4B rabbit pAb

Cat No.:ES19127

For research use only

Overview

Product Name	ANS4B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ANS4B AA range: 358-408
Specificity	This antibody detects endogenous levels of human ANS4B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Ankyrin repeat and SAM domain-containing protein 4B (Harmonin-interacting ankyrin repeat-containing protein) (Harp)
Gene Name	ANKS4B HARP
Cellular localization	Cell projection, microvillus . Localizes at the tip of microvilli (PubMed:26812018). May associate with endoplasmic reticulum membranes (By similarity). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	257629
Human Swiss-Prot Number	Q8N8V4
Alternative Names	
Background	

