



YLAT1 rabbit pAb

Cat No.:ES19131

For research use only

Overview

Product Name	YLAT1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human YLAT1 AA range: 239-289
Specificity	This antibody detects endogenous levels of human YLAT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Y+L amino acid transporter 1 (Monocyte amino acid permease 2) (MOP-2) (Solute carrier family 7 member 7) (y(+)-L-type amino acid transporter 1) (Y+LAT1) (y+LAT-1)
Gene Name	SLC7A7
Cellular localization	Basolateral cell membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9056
Human Swiss-Prot Number	Q9UM01
Alternative Names	
Background	

