



FA156 rabbit pAb

Cat No.:ES19132

For research use only

Overview

Product Name	FA156 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FA156 AA range: 6-56
Specificity	This antibody detects endogenous levels of human FA156
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein FAM156A/FAM156B (Transmembrane protein 29)
Gene Name	FAM156A TMEM29 PP12994 PRO0659; FAM156B TMEM29B
Cellular localization	Membrane ; Single-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	29057
Human Swiss-Prot Number	Q8NDB6
Alternative Names	
Background	

