



HELQ rabbit pAb

Cat No.:ES19137

For research use only

Overview

Product Name	HELQ rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HELQ AA range: 830-880
Specificity	This antibody detects endogenous levels of human HELQ
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Helicase POLQ-like (EC 3.6.4.12) (Mus308-like helicase) (POLQ-like helicase)
Gene Name	HELQ HEL308
Cellular localization	Nucleus . Chromosome . Localizes to sites of DNA damage; localizes to damaged replication forks. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	113510
Human Swiss-Prot Number	Q8TDG4
Alternative Names	
Background	

