



NUD16 rabbit pAb

Cat No.:ES19138

For research use only

Overview

Product Name	NUD16 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NUD16 AA range: 2-52
Specificity	This antibody detects endogenous levels of human NUD16
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	U8 snoRNA-decapping enzyme (EC 3.6.1.62) (Nucleoside diphosphate-linked moiety X motif 16) (Nudix motif 16) (U8 snoRNA-binding protein H29K) (m7GpppN-mRNA hydrolase)
Gene Name	NUDT16
Cellular localization	Nucleus . Nucleus, nucleoplasm . Nucleus, nucleolus . Cytoplasm . Localized predominantly in the cytoplasm (PubMed:21070968). Localized in nucleolus, and in a minor proportion in distinct foci in the nucleoplasm (By similarity). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	131870
Human Swiss-Prot Number	Q96DE0
Alternative Names	
Background	

