



H2B2C rabbit pAb

Cat No.:ES19139

For research use only

Overview

Product Name	H2B2C rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human H2B2C AA range: 82-132
Specificity	This antibody detects endogenous levels of human H2B2C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Putative histone H2B type 2-C (Histone H2B.t) (H2B/t)
Gene Name	HIST2H2BC
Cellular localization	Nucleus. Chromosome.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q6DN03
Alternative Names	
Background	

