



HRSL2 rabbit pAb

Cat No.:ES19141

For research use only

Overview

Product Name	HRSL2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HRSL2 AA range: 1-51
Specificity	This antibody detects endogenous levels of human HRSL2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	HRAS-like suppressor 2 (EC 2.3.1.-) (EC 3.1.1.-)
Gene Name	HRASLS2
Cellular localization	Cytoplasm . Membrane ; Single-pass membrane protein . Exhibits a granular pattern in the cytoplasm with preferential perinuclear localization. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	54979
Human Swiss-Prot Number	Q9NWW9
Alternative Names	
Background	

