

NUD22 rabbit pAb

Cat No.:ES19143

For research use only

Overview

Product Name	NUD22 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NUD22 AA range: 44-94
Specificity	This antibody detects endogenous levels of human NUD22
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Nucleoside diphosphate-linked moiety X motif 22 (Nudix motif 22)
Gene Name	NUDT22 PP11246
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84304
Human Swiss-Prot Number	Q9BRQ3
Alternative Names	
Background	

