



H2B1J rabbit pAb

Cat No.:ES19154

For research use only

Overview

Product Name	H2B1J rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human H2B1J AA range: 58-108
Specificity	This antibody detects endogenous levels of human H2B1J
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Histone H2B type 1-J (Histone H2B.1) (Histone H2B.r) (H2B/r)
Gene Name	HIST1H2BJ H2BFR
Cellular localization	Nucleus. Chromosome.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	8970
Human Swiss-Prot Number	P06899
Alternative Names	
Background	

