



FBX3 rabbit pAb

Cat No.:ES19156

For research use only

Overview

Product Name	FBX3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FBX3 AA range: 245-295
Specificity	This antibody detects endogenous levels of human FBX3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	F-box only protein 3
Gene Name	FBXO3 FBX3
Cellular localization	Nucleus. Colocalizes with PML at the peripheries of nuclear bodies.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	26273
Human Swiss-Prot Number	Q9UK99
Alternative Names	
Background	

