

TMOD4 rabbit pAb

Cat No.:ES19162

For research use only

Overview

Product Name	TMOD4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TMOD4 AA range: 90-140
Specificity	This antibody detects endogenous levels of human TMOD4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Tropomodulin-4 (Skeletal muscle tropomodulin) (Sk-Tmod)
Gene Name	TMOD4
Cellular localization	Cytoplasm, cytoskeleton . In myofibrils with sarcomeric structure, localizes to the pointed end of actin thin filaments (PubMed:25250574). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	29765
Human Swiss-Prot Number	Q9NZQ9
Alternative Names	
Background	

