



TEN2 rabbit pAb

Cat No.:ES19168

For research use only

Overview

Product Name	TEN2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TEN2 AA range: 2379-2429
Specificity	This antibody detects endogenous levels of human TEN2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Teneurin-2 (Ten-2) (Protein Odd Oz/ten-m homolog 2) (Tenascin-M2) (Ten-m2)
Gene Name	ODZ2 KIAA1127 TNM2
Cellular localization	Cell membrane ; Single-pass membrane protein . Endoplasmic reticulum . Golgi apparatus . Cell junction, synapse . Cell projection, dendritic spine . Cell projection, filopodium . Cell projection, growth cone . Cell junction, synapse, postsynaptic cell membrane ; Single-pass membrane protein . Cell junction, synapse, synaptosome . Colocalizes with ADGRL1 across intercellular junctions. .; [Isoform 2]: Cell membrane. Colocalizes with ADGRL1 across intercellular junctions.; [Ten-2 intracellular domain]: Nucleus, PML body .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	





ELK Biotechnology

Human Swiss-Prot Number Q9NT68
Alternative Names
Background



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C