



# RGCC rabbit pAb

Cat No.:ES19173

For research use only

## Overview

<b>Product Name</b>	RGCC rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB
<b>Species Cross-Reactivity</b>	Human;Mouse;Rat
<b>Recommended dilutions</b>	WB 1:1000-2000
<b>Immunogen</b>	Synthesized peptide derived from human RGCC AA range: 5-55
<b>Specificity</b>	This antibody detects endogenous levels of human RGCC
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C. Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	Regulator of cell cycle RGCC (Response gene to complement 32 protein) (RGC-32)
<b>Gene Name</b>	RGCC C13orf15 RGC32
<b>Cellular localization</b>	Cytoplasm. Nucleus. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome. Cytoplasmic in unstimulated cells. Nuclear after activation by complement. Associated with the centrosome during prometaphase and metaphase.
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	
<b>Human Gene ID</b>	28984
<b>Human Swiss-Prot Number</b>	Q9H4X1
<b>Alternative Names</b>	
<b>Background</b>	

