



RIM3A rabbit pAb

Cat No.:ES19175

For research use only

Overview

Product Name	RIM3A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human RIM3A AA range: 822-872
Specificity	This antibody detects endogenous levels of human RIM3A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	RIMS-binding protein 3A (RIM-BP3.A) (RIMS-binding protein 3.1) (RIM-BP3.1)
Gene Name	RIMBP3 KIAA1666 RIMBP3A
Cellular localization	Cytoplasm, cytoskeleton . In elongating spermatids, localizes to the manchette. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	85376
Human Swiss-Prot Number	Q9UFD9
Alternative Names	
Background	

