



REEP6 rabbit pAb

Cat No.:ES19180

For research use only

Overview

Product Name	REEP6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human REEP6 AA range: 110-160
Specificity	This antibody detects endogenous levels of human REEP6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Receptor expression-enhancing protein 6 (Polyposis locus protein 1-like 1)
Gene Name	REEP6 C19orf32 DP1L1
Cellular localization	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Cytoplasmic vesicle, clathrin-coated vesicle membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	92840
Human Swiss-Prot Number	Q96HR9
Alternative Names	
Background	

