



# RHBL2 rabbit pAb

Cat No.:ES19182

For research use only

## Overview

<b>Product Name</b>	RHBL2 rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB
<b>Species Cross-Reactivity</b>	Human;Mouse
<b>Recommended dilutions</b>	WB 1:1000-2000
<b>Immunogen</b>	Synthesized peptide derived from human RHBL2 AA range: 173-223
<b>Specificity</b>	This antibody detects endogenous levels of human RHBL2
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C. Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	Rhomboid-related protein 2 (RRP2) (EC 3.4.21.105) (Rhomboid-like protein 2) [Cleaved into: Rhomboid-related protein 2, N-terminal fragment (NTF); Rhomboid-related protein 2, C-terminal fragment (CTF)]
<b>Gene Name</b>	RHBDL2
<b>Cellular localization</b>	[Rhomboid-related protein 2, C-terminal fragment]: Cell membrane ; Multi-pass membrane protein .
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	
<b>Human Gene ID</b>	54933
<b>Human Swiss-Prot Number</b>	Q9NX52
<b>Alternative Names</b>	
<b>Background</b>	

