

NALP6 rabbit pAb

Cat No.:ES19189

For research use only

Overview

Product Name	NALP6 rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000 IHC-p 1:50-200;IF(paraffin section)1:50-200
Immunogen	Synthesized peptide derived from human NALP6.AA range: 20-50
Specificity	This antibody detects endogenous levels of human NALP6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	NACHT, LRR and PYD domains-containing protein 6 (Angiotensin II/vasopressin receptor) (PYRIN-containing APAF1-like protein 5)
Gene Name	NLRP6 NALP6 PYPAF5
Cellular localization	Cytoplasm, cytosol . Inflammasome . Cell membrane . Nucleus membrane .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	171389
Human Swiss-Prot Number	P59044
Alternative Names	
Background	

