



KRA56 rabbit pAb

Cat No.:ES19194

For research use only

Overview

Product Name	KRA56 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KRA56 AA range: 18-68
Specificity	This antibody detects endogenous levels of human KRA56
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Keratin-associated protein 5-6 (Keratin-associated protein 5.6) (Ultrahigh sulfur keratin-associated protein 5.6)
Gene Name	KRTAP5-6 KAP5.6 KRTAP5.6
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	440023
Human Swiss-Prot Number	Q6L8G9
Alternative Names	
Background	

