



S29P2 rabbit pAb

Cat No.:ES19197

For research use only

Overview

Product Name	S29P2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human S29P2 AA range: 105-155
Specificity	This antibody detects endogenous levels of human S29P2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Putative protein SNX29P2 (RUN domain-containing protein 2C) (Sorting nexin 29 protein pseudogene 2)
Gene Name	SNX29P2 RUNDC2C
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q8IUI4
Alternative Names	
Background	

