

HEMGN rabbit pAb

Cat No.:ES19205

For research use only

Overview

Product Name	HEMGN rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HEMGN AA range: 419-469
Specificity	This antibody detects endogenous levels of human HEMGN
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Hemogen (Erythroid differentiation-associated gene protein) (EDAG-1) (Hemopoietic gene protein) (Negative differentiation regulator protein)
Gene Name	HEMGN EDAG NDR PRO1037 PRO1620
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	55363
Human Swiss-Prot Number	Q9BXL5
Alternative Names	
Background	

