

# SEZ6 rabbit pAb

Cat No.:ES19208

For research use only

## Overview

Product Name	SEZ6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SEZ6 AA range: 108-158
Specificity	This antibody detects endogenous levels of human SEZ6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Seizure protein 6 homolog (SEZ-6) (hSEZ-6)
Gene Name	SEZ6
Cellular localization	Cell membrane ; Single-pass type I membrane protein . Localized on dendrites and in the synaptic and postsynaptic fraction. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	124925
Human Swiss-Prot Number	Q53EL9
Alternative Names	
Background	

