



NDK5 rabbit pAb

Cat No.:ES19209

For research use only

Overview

Product Name	NDK5 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NDK5 AA range: 25-75
Specificity	This antibody detects endogenous levels of human NDK5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Nucleoside diphosphate kinase homolog 5 (NDK-H 5) (NDP kinase homolog 5) (Inhibitor of p53-induced apoptosis-beta) (IPIA-beta) (Testis-specific nm23 homolog) (nm23-H5)
Gene Name	NME5
Cellular localization	Cell projection, cilium .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	8382
Human Swiss-Prot Number	P56597
Alternative Names	
Background	

