



TSN5 rabbit pAb

Cat No.:ES19210

For research use only

Overview

Product Name	TSN5 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TSN5 AA range: 113-163
Specificity	This antibody detects endogenous levels of human TSN5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Tetraspanin-5 (Tspan-5) (Tetraspan NET-4) (Transmembrane 4 superfamily member 9)
Gene Name	TSPAN5 TM4SF9
Cellular localization	Cell membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	10098
Human Swiss-Prot Number	P62079
Alternative Names	
Background	

