



RBG1L rabbit pAb

Cat No.:ES19214

For research use only

Overview

Product Name	RBG1L rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human RBG1L AA range: 639-689
Specificity	This antibody detects endogenous levels of human RBG1L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Rab GTPase-activating protein 1-like
Gene Name	RABGAP1L HHL KIAA0471
Cellular localization	Cytoplasmic vesicle . Early endosome . Golgi apparatus . Colocalizes on endosomes partially with EEA1 (PubMed:16923123). Colocalizes and cotransports on motile vesicles with ANK2 (By similarity). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9910
Human Swiss-Prot Number	Q5R372
Alternative Names	
Background	

