



F159A rabbit pAb

Cat No.:ES19216

For research use only

Overview

Product Name	F159A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human F159A AA range: 88-138
Specificity	This antibody detects endogenous levels of human F159A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Membrane protein FAM159A
Gene Name	FAM159A UNQ2783/PRO7171
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	348378
Human Swiss-Prot Number	Q6UWV7
Alternative Names	
Background	

