

# UN13C rabbit pAb

Cat No.:ES19217

For research use only

## Overview

Product Name	UN13C rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human UN13C AA range: 1996-2046
Specificity	This antibody detects endogenous levels of human UN13C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein unc-13 homolog C (Munc13-3)
Gene Name	UNC13C
Cellular localization	Cytoplasm . Membrane ; Peripheral membrane protein . Cell junction, synapse, presynaptic cell membrane ; Peripheral membrane protein . Localized to presynaptic structures. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	440279
Human Swiss-Prot Number	Q8NB66
Alternative Names	
Background	

