



VKIND rabbit pAb

Cat No.:ES19219

For research use only

Overview

Product Name	VKIND rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human VKIND AA range: 1464-1514
Specificity	This antibody detects endogenous levels of human VKIND
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein very KIND (Cerebral protein 9) (KIND domain-containing protein 1) (Kinase non-catalytic C-lobe domain-containing protein 1) (Ras-GEF domain-containing family member 2)
Gene Name	KNDC1 C10orf23 KIAA1768 RASGEF2 VKIND hucep-9
Cellular localization	Cell projection, dendrite . Perikaryon .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	85442
Human Swiss-Prot Number	Q76NI1
Alternative Names	
Background	

