

BTBD8 rabbit pAb

Cat No.:ES19224

For research use only

Overview

Product Name	BTBD8 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human BTBD8 AA range: 224-274
Specificity	This antibody detects endogenous levels of human BTBD8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	BTB/POZ domain-containing protein 8
Gene Name	BTBD8
Cellular localization	Nucleus . Localized to nucleus in fetal cells.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	284697
Human Swiss-Prot Number	Q5XKL5
Alternative Names	
Background	

