

TRNK1 rabbit pAb

Cat No.:ES19228

For research use only

Overview

Product Name	TRNK1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TRNK1 AA range: 2101-2151
Specificity	This antibody detects endogenous levels of human TRNK1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	TPR and ankyrin repeat-containing protein 1 (Lupus brain antigen 1 homolog)
Gene Name	TRANK1 KIAA0342 LBA1
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9881
Human Swiss-Prot Number	O15050
Alternative Names	
Background	

