

# CQ10A rabbit pAb

Cat No.:ES19232

For research use only

## Overview

Product Name	CQ10A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CQ10A AA range: 136-186
Specificity	This antibody detects endogenous levels of human CQ10A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Coenzyme Q-binding protein COQ10 homolog A, mitochondrial
Gene Name	COQ10A UNQ6192/PRO20219
Cellular localization	Mitochondrion inner membrane ; Peripheral membrane protein ; Matrix side .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	93058
Human Swiss-Prot Number	Q96MF6
Alternative Names	
Background	

