



FA83F rabbit pAb

Cat No.:ES19236

For research use only

Overview

Product Name	FA83F rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FA83F AA range: 324-374
Specificity	This antibody detects endogenous levels of human FA83F
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Protein FAM83F
Gene Name	FAM83F
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	113828
Human Swiss-Prot Number	Q8NEG4
Alternative Names	
Background	

