

# VP13B rabbit pAb

Cat No.:ES19238

For research use only

## Overview

Product Name	VP13B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human VP13B AA range: 3593-3643
Specificity	This antibody detects endogenous levels of human VP13B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Vacuolar protein sorting-associated protein 13B (Cohen syndrome protein 1)
Gene Name	VPS13B CHS1 COH1 KIAA0532
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	157680
Human Swiss-Prot Number	Q7Z7G8
Alternative Names	
Background	

