

JAGN1 rabbit pAb

Cat No.:ES19244

For research use only

Overview

Product Name	JAGN1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human JAGN1 AA range: 111-161
Specificity	This antibody detects endogenous levels of human JAGN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein jagunal homolog 1
Gene Name	JAGN1
Cellular localization	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84522
Human Swiss-Prot Number	Q8N5M9
Alternative Names	
Background	

