



CNTP5 rabbit pAb

Cat No.:ES19245

For research use only

Overview

Product Name	CNTP5 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CNTP5 AA range: 765-815
Specificity	This antibody detects endogenous levels of human CNTP5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Contactin-associated protein-like 5 (Cell recognition molecule Caspr5)
Gene Name	CNTNAP5 CASPR5
Cellular localization	Membrane ; Single-pass type I membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	129684
Human Swiss-Prot Number	Q8WYK1
Alternative Names	
Background	

