



S4A8 rabbit pAb

Cat No.:ES19248

For research use only

Overview

Product Name	S4A8 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human S4A8 AA range: 1012-1062
Specificity	This antibody detects endogenous levels of human S4A8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Electroneutral sodium bicarbonate exchanger 1 (Electroneutral Na(+)-driven Cl-HCO ₃ exchanger) (Solute carrier family 4 member 8) (k-NBC3)
Gene Name	SLC4A8 KIAA0739 NBC NBC3 NDCBE1
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9498
Human Swiss-Prot Number	Q2Y0W8
Alternative Names	
Background	

