

ASXL1 rabbit pAb

Cat No.: ES19255

For research use only

Overview

Product Name ASXL1 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human; Mouse Recommended dilutions WB 1:1000-2000

Immunogen Synthesized peptide derived from human ASXL1 AA

range: 431-481

Specificity This antibody detects endogenous levels of human

ASXL1

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Putative Polycomb group protein ASXL1 (Additional

sex combs-like protein 1)

Gene Name ASXL1 KIAA0978

Cellular localization Nucleus.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/ml

Observed band

Human Gene ID 171023 Human Swiss-Prot Number Q8IXJ9

Alternative Names

Background

