



LRC45 rabbit pAb

Cat No.:ES19256

For research use only

Overview

Product Name	LRC45 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LRC45 AA range: 483-533
Specificity	This antibody detects endogenous levels of human LRC45
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Leucine-rich repeat-containing protein 45
Gene Name	LRRC45
Cellular localization	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Localizes to the proteinaceous linker between the proximal ends of the centrioles. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	201255
Human Swiss-Prot Number	Q96CN5
Alternative Names	
Background	

