

# RN152 rabbit pAb

Cat No.:ES19257

For research use only

## Overview

Product Name	RN152 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human RN152 AA range: 150-200
Specificity	This antibody detects endogenous levels of human RN152
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	E3 ubiquitin-protein ligase RNF152 (EC 6.3.2.-) (RING finger protein 152)
Gene Name	RNF152
Cellular localization	Lysosome membrane ; Single-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	220441
Human Swiss-Prot Number	Q8N8N0
Alternative Names	
Background	

