



# S26A1 rabbit pAb

Cat No.:ES19259

For research use only

## Overview

Product Name	S26A1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human S26A1 AA range: 502-552
Specificity	This antibody detects endogenous levels of human S26A1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Sulfate anion transporter 1 (SAT-1) (Solute carrier family 26 member 1)
Gene Name	SLC26A1 SAT1
Cellular localization	Cell membrane ; Multi-pass membrane protein . Basolateral cell membrane .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	10861
Human Swiss-Prot Number	Q9H2B4
Alternative Names	
Background	

