

CG061 rabbit pAb

Cat No.:ES19260

For research use only

Overview

Product Name	CG061 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CG061 AA range: 45-95
Specificity	This antibody detects endogenous levels of human CG061
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Uncharacterized protein C7orf61
Gene Name	C7orf61
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	402573
Human Swiss-Prot Number	Q8IZ16
Alternative Names	
Background	

