



DEDD2 rabbit pAb

Cat No.:ES19264

For research use only

Overview

Product Name	DEDD2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DEDD2 AA range: 226-276
Specificity	This antibody detects endogenous levels of human DEDD2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	DNA-binding death effector domain-containing protein 2 (DED-containing protein FLAME-3) (FADD-like anti-apoptotic molecule 3)
Gene Name	DEDD2 FLAME3 PSEC0004
Cellular localization	Nucleus, nucleolus . Nuclear, accumulated in subnuclear structures resembling nucleoli.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	162989
Human Swiss-Prot Number	Q8WXF8
Alternative Names	
Background	

