

# KR103 rabbit pAb

Cat No.:ES19265

For research use only

## Overview

Product Name	KR103 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KR103 AA range: 27-77
Specificity	This antibody detects endogenous levels of human KR103
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Keratin-associated protein 10-3 (High sulfur keratin-associated protein 10.3) (Keratin-associated protein 10.3) (Keratin-associated protein 18-3) (Keratin-associated protein 18.3)
Gene Name	KRTAP10-3 KAP10.3 KAP18-3 KRTAP10.3 KRTAP18-3 KRTAP18.3
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	386682
Human Swiss-Prot Number	P60369
Alternative Names	
Background	

