



DCAF8 rabbit pAb

Cat No.:ES19271

For research use only

Overview

Product Name	DCAF8 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DCAF8 AA range: 65-115
Specificity	This antibody detects endogenous levels of human DCAF8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	DDB1- and CUL4-associated factor 8 (WD repeat-containing protein 42A)
Gene Name	DCAF8 H326 WDR42A
Cellular localization	Nucleus . Cytoplasm . It shuttles between the nucleus and the cytoplasm. Nuclear import is mediated by KPNA1 and KPNB1 under the regulation of nuclear GTPase RAN. Nuclear export to the cytoplasm is XPO1 dependent. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	50717
Human Swiss-Prot Number	Q5TAQ9
Alternative Names	
Background	

